

File Name: 5_D_SB_PRIME.LHC
File Location: C:\ProgramData\BioTek\Liquid Handling Control
2.22\Protocols\automation_v3.1
Last Saved: 31/05/2021 14:07:23

LHC Version: 2.19.1

Instrument: MultiFloFX
Port: USB MultiFloFX sn:19041612
Settings: Non-Autoclavable Syringes
16-tube 7° Syringe manifold
Peri-pump Dispenser
Peri-pump Dispenser 2
Random Access Dispenser

Plate Type: 384 Well Plate
Protocol Name: P2_RAD-384
Protocol Version: 1.0
Archive Revision: 46
Comments: BioTek provided protocol to demonstrate dispensing into a 384 well plate using
RAD technology on instruments with an integrated peripump.

Step Details: OBS! catch the PFA waste in a container.
S-Prime
Syringe: B
Volume: 5000 µL
Flow Rate: 5
Number of Cycles: 4
Pump Delay: 0 msec
Submerge tips: No

<end>