

File Name: 2_D_P1_40ul_mito.LHC
File Location: C:\ProgramData\BioTek\Liquid Handling Control
2.22\Protocols\automation_v3.1
Last Saved: 07/10/2021 09:28:28

LHC Version: 2.19.1

Instrument: MultiFloFX
Port: USB MultiFloFX sn:19041612
Settings: Non-Autoclavable Syringes
16-tube 7° Syringe manifold
Peri-pump Dispenser
Peri-pump Dispenser 2
Random Access Dispenser

Plate Type: 384 Well Plate
Protocol Name: P1_RAD-384
Protocol Version: 1.0
Archive Revision: 75
Comments: OBS!! when 5ul cassette attached: 160 ul= 80 ul

no predispensing here, the prime part of the other mito protocol does the
predispensing
Step Details: P-Dispense
Peri-pump: Primary
Volume: 40 µL/tube
Flow Rate: Low
Required Cassette: Any
Z Offset: 333 steps (15.22 mm above carrier)
X Offset: -20 steps (0.91 mm left of center)
Y Offset: 0 steps (center of well)
Pre-dispense: No
Columns: 11111111111111111111111111111111
Rows: 11

<end>