

File Name: 2\_D\_P1\_purge\_then\_prime.LHC  
File Location: C:\ProgramData\BioTek\Liquid Handling Control  
2.22\Protocols\automation\_v3.1  
Last Saved: 07/10/2021 09:28:26

LHC Version: 2.19.1

Instrument: MultiFloFX  
Port: USB MultiFloFX sn:19041612  
Settings: Non-Autoclavable Syringes  
16-tube 7° Syringe manifold  
Peri-pump Dispenser  
Peri-pump Dispenser 2  
Random Access Dispenser

Plate Type: 384 Well Plate  
Protocol Name: P1\_RAD-384  
Protocol Version: 1.0  
Archive Revision: 75  
Comments: OBS!! when 5ul cassette attached: 160 ul= 80 ul

Step Details: priming = predispensing  
P-Purge  
Peri-pump: Primary  
Volume: 2600 µL/tube  
Flow Rate: High  
Required Cassette: Any

P-Prime  
Peri-pump: Primary  
Volume: 1400 µL/tube  
Flow Rate: Low  
Required Cassette: Any

<end>