

File Name: 3\_D\_SA\_PRIME.LHC  
File Location: C:\ProgramData\BioTek\Liquid Handling Control  
2.22\Protocols\automation\_v3.1  
Last Saved: 31/05/2021 14:07:13

LHC Version: 2.19.1

Instrument: MultiFloFX  
Port: USB MultiFloFX sn:19041612  
Settings: Non-Autoclavable Syringes  
16-tube 7° Syringe manifold  
Peri-pump Dispenser  
Peri-pump Dispenser 2  
Random Access Dispenser

Plate Type: 384 Well Plate  
Protocol Name: P2\_RAD-384  
Protocol Version: 1.0  
Archive Revision: 45  
Comments: BioTek provided protocol to demonstrate dispensing into a 384 well plate using  
RAD technology on instruments with an integrated peripump.

Step Details: OBS! catch the PFA waste in a container.  
S-Prime  
Syringe: A  
Volume: 5000 µL  
Flow Rate: 5  
Number of Cycles: 4  
Pump Delay: 0 msec  
Submerge tips: No

<end>