

File Name: 8_D_P2_20ul_stains.LHC
File Location: C:\ProgramData\BioTek\Liquid Handling Control
2.22\Protocols\automation_v3.1
Last Saved: 07/10/2021 09:31:42

LHC Version: 2.19.1

Instrument: MultiFloFX
Port: USB MultiFloFX sn:19041612
Settings: Non-Autoclavable Syringes
16-tube 7° Syringe manifold
Peri-pump Dispenser
Peri-pump Dispenser 2
Random Access Dispenser

Plate Type: 384 Well Plate
Protocol Name: P2_RAD-384
Protocol Version: 1.0
Archive Revision: 52
Comments: 20 = 10 ul

Step Details: no predispensing here, priming is done in the other protocol
P-Dispense
Peri-pump: Secondary
Volume: 20 µL/tube
Flow Rate: Medium
Required Cassette: Any
Z Offset: 333 steps (15.22 mm above carrier)
X Offset: 0 steps (center of well)
Y Offset: 0 steps (center of well)
Pre-dispense: No
Columns: 111111111111111111111111
Rows: 11

<end>