

File Name: 4\_D\_SA\_384\_80ul\_PFA.LHC  
File Location: C:\ProgramData\BioTek\Liquid Handling Control  
2.22\Protocols\automation\_v3.1  
Last Saved: 01/09/2021 13:27:47

LHC Version: 2.19.1

Instrument: MultiFloFX  
Port: USB MultiFloFX sn:19041612  
Settings: Non-Autoclavable Syringes  
16-tube 7° Syringe manifold  
Peri-pump Dispenser  
Peri-pump Dispenser 2  
Random Access Dispenser

Plate Type: 384 Well Plate  
Protocol Name: P2\_RAD-384  
Protocol Version: 1.0  
Archive Revision: 44  
Comments: BioTek provided protocol to demonstrate dispensing into a 384 well plate using  
RAD technology on instruments with an integrated peripump.

Step Details: OBS! catch the PFA waste in a container.  
S-Dispense  
Syringe: A  
Volume: 80 µL/well  
Flow Rate: 1  
Z Offset: 333 steps (15.22 mm above carrier)  
X Offset: 0 steps (center of well)  
Y Offset: 0 steps (center of well)  
Pre-dispense: No  
Pump Delay: 0 msec  
Columns: 111111111111111111111111

<end>