

File Name: 8_D_P2_purge_then_prime.LHC
File Location: C:\ProgramData\BioTek\Liquid Handling Control
2.22\Protocols\automation_v3.1
Last Saved: 07/10/2021 09:31:12

LHC Version: 2.19.1

Instrument: MultiFloFX
Port: USB MultiFloFX sn:19041612
Settings: Non-Autoclavable Syringes
16-tube 7° Syringe manifold
Peri-pump Dispenser
Peri-pump Dispenser 2
Random Access Dispenser

Plate Type: 384 Well Plate
Protocol Name: P2_RAD-384
Protocol Version: 1.0
Archive Revision: 52
Comments: 20 = 10 ul

Step Details: priming part does the predispensing
P-Purge
Peri-pump: Secondary
Volume: 1500 µL/tube
Flow Rate: High
Required Cassette: Any

P-Prime
Peri-pump: Secondary
Volume: 460 µL/tube
Flow Rate: Medium
Required Cassette: Any

<end>